

CAB0136

Product Information

Product SKU:	CAB0136	Gene ID:	27037	Size:	20uL, 100uL
Clone No:	ARC2498	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	75kDa	Conjugate:	Unmodified
Calculated MW:	69kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is of unknown function. However, it is orthologous to the mouse Trmt2a gene and contains an RNA methyltransferase domain. Expression of this gene varies during the cell cycle, with aberrant expression being a possible biomarker in certain breast cancers. Several transcript variants encoding two different isoforms have been found for this gene.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: HTF9C; TRMT2A/HTF9C

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TRMT2A/HTF9C (Q8IZ69).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.