

CAB0341

Product Information

Product SKU:	CAB0341	Gene ID:	23016	Size:	20uL, 100uL
Clone No:	ARC2501	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	37kDa	Conjugate:	Unconjugated
Calculated MW:	32kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable 3'-5'-exoribonuclease activity and RNA binding activity. Predicted to be involved in RNA metabolic process. Part of exosome (RNase complex).

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: p8; EAP1; RRP42; Rrp42p; hRrp42p; EXOSC7

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 192-291 of human EXOSC7 (Q15024).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.