

CAB0381

Product Information

Product SKU: CAB0381	Gene ID: 891	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: 60kDa	Conjugate: Unconjugated
Calculated MW: 48kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a regulatory protein involved in mitosis. The gene product complexes with p34(cdc2) to form the maturation-promoting factor (MPF). The encoded protein is necessary for proper control of the G2/M transition phase of the cell cycle.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: CCNB; Cyclin B1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 350-433 of human Cyclin B1 (NP_114172.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.