

## CAB0455

---

### Product Information

<b>Product SKU:</b>	CAB0455	<b>Gene ID:</b>	6317	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2508	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	Refer to figures	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	45kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables cysteine-type endopeptidase inhibitor activity; protease binding activity; and virus receptor activity. Involved in several processes, including autocrine signaling; paracrine signaling; and regulation of cellular protein metabolic process. Located in cytoplasmic vesicle; extracellular exosome; and nucleus.

**Recommended Dilution:** IHC-P,1:50 - 1:200

**Synonyms:** SCC; T4-A; SCCA1; SSSA1; SCCA-1; HsT1196; SCCA-PD; SSSA1/SerpinB3

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 301-390 of human SSSA1/SerpinB3 (P29508).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.