

CAB0501

Product Information

Product SKU:	CAB0501	Gene ID:	23016	Size:	20uL, 100uL
Clone No:	ARC2510	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	37kDa	Conjugate:	Unmodified
Calculated MW:	32kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable 3'-5'-exoribonuclease activity and RNA binding activity. Predicted to be involved in RNA metabolic process. Part of exosome (RNase complex).

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: p8; EAP1; RRP42; Rrp42p; hRrp42p; EXOSC7

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human EXOSC7 (Q15024).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.