

CAB0621

Product Information

Product SKU:	CAB0621	Gene ID:	23358	Size:	20uL, 100uL
Clone No:	ARC2526	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	295kDa	Conjugate:	Unmodified
Calculated MW:	294kDa	Isotype:	IgG

Immunogen Information

Background: Modification of cellular proteins by ubiquitin is an essential regulatory mechanism controlled by the coordinated action of multiple ubiquitin-conjugating and deubiquitinating enzymes. USP24 belongs to a large family of cysteine proteases that function as deubiquitinating enzymes (Quesada et al., 2004 [PubMed 14715245]).

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: USP24; ubiquitin specific peptidase 24

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 2500-2600 of human USP24 (Q9UPU5).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.