

CAB0623

Product Information

Product SKU:	CAB0623	Gene ID:	29079	Size:	20uL, 100uL
Clone No:	ARC2528	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	36kDa	Conjugate:	Unmodified
Calculated MW:	30kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a component of the Mediator complex. The Mediator complex interacts with DNA-binding gene-specific transcription factors to modulate transcription by RNA polymerase II. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: ARC36; VDRIP; DRIP36; TRAP36; HSPC126; MED4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 150-250 of human MED4 (Q9NPJ6).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.