

CAB0664

Product Information

Product SKU:	CAB0664	Gene ID:	10449	Size:	20uL, 100uL
Clone No:	ARC2533	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	Unmodified
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background: The encoded protein catalyzes the last step of the mitochondrial fatty acid beta-oxidation spiral. Unlike most mitochondrial matrix proteins, it contains a non-cleavable amino-terminal targeting signal.

Recommended Dilution: WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

Synonyms: DSAEC; ACAA2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 298-397 of human ACAA2 (P42765).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.