

CAB0723

Product Information

Product SKU:	CAB0723	Gene ID:	4100	Size:	20uL, 100uL
Clone No:	ARC1839	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	51kDa	Conjugate:	Unconjugated
Calculated MW:	34kDa	Isotype:	IgG

Immunogen Information

Background: This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: CT1.1; MAGE1; MAGE-1/MAGE1A

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 200-309 of human MAGE-1/MAGE1A (P43355).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.