

CAB0767

Product Information

Product SKU: CAB0767	Gene ID: 301	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Rat

Additional Information

Observed MW: 34kDa	Conjugate: Unconjugated
Calculated MW: 39kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a membrane-localized protein that binds phospholipids. This protein inhibits phospholipase A2 and has anti-inflammatory activity. Loss of function or expression of this gene has been detected in multiple tumors.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: ANX1; LPC1; ANXA1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-346 of human ANXA1 (NP_000691.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.