

CAB0813

Product Information

Product SKU:	CAB0813	Gene ID:	4597	Size:	20uL, 100uL
Clone No:	ARC2541	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	43kDa	Conjugate:	Unmodified
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: The enzyme mevalonate pyrophosphate decarboxylase catalyzes the conversion of mevalonate pyrophosphate into isopentenyl pyrophosphate in one of the early steps in cholesterol biosynthesis. It decarboxylates and dehydrates its substrate while hydrolyzing ATP.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: MPD; MDDase; POROK7; FP17780; MVD

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 250-350 of human MVD (P53602).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.