

CAB0914

Product Information

Product SKU:	CAB0914	Gene ID:	56998	Size:	20uL, 100uL
Clone No:	ARC2551	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	11kDa	Conjugate:	Unmodified
Calculated MW:	9kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene binds CTNNB1 and prevents interaction between CTNNB1 and TCF family members. The encoded protein is a negative regulator of the Wnt signaling pathway. Two transcript variants encoding the same protein have been found for this gene.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: ICAT; CTNNBIP1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-81 of human CTNNBIP1 (Q9NSA3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.