

CAB10341

Product Information

Product SKU:	CAB10341	Gene ID:	729440	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	57kDa	Conjugate:	Unconjugated
Calculated MW:	57kDa	Isotype:	IgG

Immunogen Information

Background: Enables identical protein binding activity and microtubule binding activity. Involved in centriole assembly and mitotic spindle assembly. Located in centriolar satellite and ciliary basal body. Part of centriolar subdistal appendage.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: VFL3; hVFL3; CCDC61

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 352-531 of human CCDC61 (NP_001073871.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.