

MEIOC Rabbit Polyclonal Antibody



CAB10522

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

107kDa

Calculated MW:

102kDa/108kDa

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

268491

Uniprot

Synonyms:

MEIOC; C17orf104; Gm1564

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of mouse MEIOC (NP_001121048.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
