

# VGLL4 Rabbit Polyclonal Antibody



CAB10578

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

31kDa

### Calculated MW:

30kDa

### Applications:

WB

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

232334

#### Uniprot

#### Synonyms:

VGL-4; BC048841; VGLL4

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-287 of mouse VGLL4 (NP\_808351.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---