

## CAB11861

---

### Product Information

<b>Product SKU:</b>	CAB11861	<b>Gene ID:</b>	10763	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	150kDa/177kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	177kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells.

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

**Synonyms:** Nbla00170; Nestin

**Purification Method:** Affinity purification

**Immunogen:** Recombinant Protein corresponding to a sequence within amino acids 1392-1621 of human NES (NP\_006608.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.