

C5orf52 Rabbit Polyclonal Antibody



CAB11987

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

18kDa

Calculated MW:

10kDa/19kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

67430

Uniprot

Q9CR34

Synonyms:

C5orf52; 4921536K21Rik; AV043724


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 20-168 of mouse C5orf52 (NP_080426.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - C5orf52 Rabbit pAb (CAB11987)