

CAB13607

Product Information

Product SKU:	CAB13607	Gene ID:	245711	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background: Enables protein kinase activator activity and protein kinase binding activity. Involved in several processes, including G1/S transition of mitotic cell cycle; positive regulation of cell population proliferation; and positive regulation of protein kinase activity. Located in nucleoplasm.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: SPY1; SPDY1; RINGO3; RINGOA; SPDYA

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human SPDYA (NP_877433.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.