

SKIV2L2 Rabbit Polyclonal Antibody



CAB13921

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

117kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23517

Uniprot

P42285

Synonyms:

SKIV2L2; Dob1; KIAA0052; Mtr4; fSAP118

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-445 of human SKIV2L2 (NP_056175.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
