

C6orf150 Rabbit Polyclonal Antibody



CAB14228

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

59kDa

Calculated MW:

49kDa/58kDa

Applications:

WB IHC

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

115004

Uniprot

Q8N884

Synonyms:

MB21D1; C6orf150; cGAS; h-cGAS; CGAS

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human C6orf150 (NP_612450.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
