

CAB14230

Product Information

Product SKU: CAB14230	Gene ID: 115209	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 37kDa	Conjugate: Unconjugated
Calculated MW: 60kDa	Isotype: IgG

Immunogen Information

Background: Enables metalloendopeptidase activity. Involved in several processes, including HRI-mediated signaling; proteolysis; and regulation of mitochondrion organization. Located in mitochondrial inner membrane.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: DAB1; MPRP1; MPRP-1; YKR087C; ZMPOMA1; peptidase; 2010001O09Rik; OMA1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 360-524 of human OMA1 (NP_660286.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.