

CAB14342

Product Information

Product SKU: CAB14342	Gene ID: 64077	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 29kDa	Conjugate: Unconjugated
Calculated MW: 29kDa	Isotype: IgG

Immunogen Information

Background: Enables inorganic diphosphatase activity and protein homodimerization activity. Involved in phosphate-containing compound metabolic process. Located in cytosol and nuclear speck.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: HDHD2B; LHPP

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human LHPP (NP_071409.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.