

## CAB1510

---

### Product Information

<b>Product SKU:</b>	CAB1510	<b>Gene ID:</b>	10081	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2564	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	55kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	55kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a 59 kDa protein that is associated with the U11 small nuclear ribonucleoprotein (snRNP), which is a component of the minor U12-type spliceosome responsible for catalyzing pre-mRNA splicing of U12-type introns.

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

**Synonyms:** 59K; ES18; HES18; PDCD7

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PDCD7 (Q8N8D1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.