

CAB1534

Product Information

Product SKU:	CAB1534	Gene ID:	10342	Size:	20uL, 100uL
Clone No:	ARC1882	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	55kDa	Conjugate:	Unconjugated
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: There are several documented fusion oncoproteins encoded partially by this gene. This gene also participates in several oncogenic rearrangements resulting in anaplastic lymphoma and mixoid chondrosarcoma, and may play a role in the NF-kappaB pathway. Multiple transcript variants have been found for this gene.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: TF6; HMSNP; SPG57; TRKT3; TRK fused gene (TFG)

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 301-400 of human TRK fused gene (TFG) (Q92734).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.