

RSG1 Rabbit Polyclonal Antibody



CAB15510

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

28kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

79363

Uniprot

Q9BU20

Synonyms:

RSG1; C1orf89

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human RSG1 (NP_112169.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
