

## CAB15816

---

### Product Information

<b>Product SKU:</b>	CAB15816	<b>Gene ID:</b>	26958	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	105kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	98kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable structural molecule activity. Involved in intra-Golgi vesicle-mediated transport. Part of COPI vesicle coat.

**Recommended Dilution:** WB,1:500 - 1:2000 IHC-P,1:50 - 1:100 IF/ICC,1:50 - 1:100

**Synonyms:** 2-COP; gamma-2-COP; COPG2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 542-871 of human COPG2 (NP\_036265.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.