

CAB15900

Product Information

Product SKU:	CAB15900	Gene ID:	80777	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	20kDa	Conjugate:	Unconjugated
Calculated MW:	17kDa	Isotype:	IgG

Immunogen Information

Background: Enables heme binding activity. Contributes to nitrite reductase (NO-forming) activity. Involved in nitric oxide biosynthetic process. Located in membrane. Part of nitric-oxide synthase complex.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: OMB5; CYB5-M; CYPB5M; CYB5B

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 16-150 of human CYB5B (NP_085056.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.