

# COL14A1 Rabbit Polyclonal Antibody



**CAB16443**

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

### Calculated MW:

183kDa/191kDa/193kDa

### Applications:

IHC IF

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

IHC 1:50 - 1:200 IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

This gene encodes the alpha chain of type XIV collagen, a member of the FACIT (fibril-associated collagens with interrupted triple helices) collagen family. Type XIV collagen interacts with the fibril surface and is involved in the regulation of fibrillogenesis.

## Immunogen information

### Gene ID:

7373

### Uniprot

Q05707

### Synonyms:

COL14A1; UND

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1060-1360 of human COL14A1 (NP\_066933.1).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---