

# FASTKD2 Rabbit Polyclonal Antibody



CAB16499

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

74kDa/81kDa

### Applications:

WB

### Reactivity:

Human

## Protein Background

This gene encodes a protein that is localized in the mitochondrial inner compartment and that may play a role in mitochondrial apoptosis. Nonsense mutations have been reported to result in cytochrome c oxidase deficiency.

## Immunogen information

### Gene ID:

22868

### Uniprot

Q9NYY8

### Synonyms:

FASTKD2; KIAA0971

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 439-710 of human FASTKD2 (NP\_055744.2).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---