

## CAB16539

---

### Product Information

<b>Product SKU:</b> CAB16539	<b>Gene ID:</b> 54819	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Mouse

---

### Additional Information

<b>Observed MW:</b> 21kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 21kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Predicted to enable nucleic acid binding activity and zinc ion binding activity.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** ZCCHC10

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human ZCCHC10 (NP\_001287745.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.