

H2A Rabbit Polyclonal Antibody



CAB16671

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

15kDa

Applications:

WB

Reactivity:

Oryza sativa

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4334071

Uniprot

Q84MP7

Synonyms:

OJ1365_D05.1; H2A

Immunogen:

A synthetic peptide of Oryza sativa H2A.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
