

CAB17158

Product Information

Product SKU:	CAB17158	Gene ID:	51621	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	45kDa	Conjugate:	-
Calculated MW:	31kDa	Isotype:	IgG

Immunogen Information

Background: KLF13 belongs to a family of transcription factors that contain 3 classical zinc finger DNA-binding domains consisting of a zinc atom tetrahedrally coordinated by 2 cysteines and 2 histidines (C2H2 motif). These transcription factors bind to GC-rich sequences and related GT and CACCC boxes (Schoy et al., 2000 [PubMed 11087666]).

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: BTEB3; FKLF2; NSLP1; RFLAT1; RFLAT-1; KLF13

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human KLF13 (NP_057079.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.