

## CAB17348

---

### Product Information

<b>Product SKU:</b>	CAB17348	<b>Gene ID:</b>	29801	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

### Additional Information

<b>Observed MW:</b>	81kDa	<b>Conjugate:</b>	-
<b>Calculated MW:</b>	81kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a four transmembrane protein that is a member of the zinc finger DHHC domain-containing protein family. The encoded protein may function as a palmitoyltransferase. Defects in this gene may be associated with a susceptibility to schizophrenia. Alternate splicing of this gene results in multiple transcript variants. A pseudogene of this gene is found on chromosome 22.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** DHHC8; ZNF378; ZDHHCL1; ZDHHC8

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 240-380 of human ZDHHC8 (NP\_037505.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.