

## CAB17556

---

### Product Information

<b>Product SKU:</b>	CAB17556	<b>Gene ID:</b>	4345998	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Oryza sativa

---

### Additional Information

<b>Observed MW:</b>	52kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	42kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Trans-acting factor that binds specifically to the consensus nucleotide sequence 5'-TNCGTACAA-3' (By similarity). May be involved in panicle development.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** WFP; IPA1; SPL14; OsSPL14

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Oryza sativa IPA1.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.