

# VPRBP Rabbit Polyclonal Antibody



CAB17596

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

169kDa

### Calculated MW:

118kDa/168kDa/169kDa

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

9730

#### Uniprot

Q9Y4B6

#### Synonyms:

DCAF1; RIP; VPRBP

#### Immunogen:


A synthetic peptide of Human VPRBP.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - VPRBP Polyclonal Antibody (CAB17596)