

BMPR1B Rabbit Polyclonal Antibody



CAB17858

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

80kDa

Calculated MW:

56kDa/60kDa

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a member of the bone morphogenetic protein (BMP) receptor family of transmembrane serine/threonine kinases. The ligands of this receptor are BMPs, which are members of the TGF-beta superfamily. BMPs are involved in endochondral bone formation and embryogenesis. These proteins transduce their signals through the formation of heteromeric complexes of 2 different types of serine (threonine) kinase receptors: type I receptors of about 50-55 kD and type II receptors of about 70-80 kD. Type II receptors bind ligands in the absence of type I receptors, but they require their respective type I receptors for signaling, whereas type I receptors require their respective type II receptors for ligand binding. Mutations in this gene have been associated with primary pulmonary hypertension. Several transcript variants encoding two different isoforms have been found for this gene.

Immunogen information

Gene ID:

658

Uniprot

O00238

Synonyms:

ALK-6; ALK6; AMDD; BDA1D; BDA2; CDw293; BMPR1B


Immunogen:

Recombinant protein of human BMPR1B.


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images




Western blot analysis of extracts of Rat brain, using BMPR1B Rabbit pAb (CAB17858) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (CABM00021). Exposure time: 3min.



Immunohistochemistry of paraffin-embedded human mammary cancer using BMPR1B antibody (CAB17858) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human placenta using BMPR1B antibody (CAB17858) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded rat ovary using BMPR1B antibody (CAB17858) at dilution of 1:200 (40x lens).