

CAB17980

Product Information

Product SKU:	CAB17980	Gene ID:	55146	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	-
Calculated MW:	40kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable protein-cysteine S-palmitoyltransferase activity. Predicted to be involved in peptidyl-L-cysteine S-palmitoylation and protein targeting to membrane. Located in Golgi apparatus and endoplasmic reticulum.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: ZNF374; ZDHHC4

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 110-210 of human ZDHHC4 (NP_060576.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.