

CAB18114

Product Information

Product SKU:	CAB18114	Gene ID:	25921	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	78kDa	Conjugate:	-
Calculated MW:	78kDa	Isotype:	IgG

Immunogen Information

Background: Enables palmitoyltransferase activity. Involved in positive regulation of pattern recognition receptor signaling pathway and positive regulation of protein localization to plasma membrane. Acts upstream of or within protein palmitoylation. Located in phagocytic vesicle and plasma membrane.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: DHHC5; ZNF375; ZDHHC5

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 130-230 of human ZDHHC5 (NP_056272.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.