

CAB18179

Product Information

Product SKU: CAB18179	Gene ID: 8464	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: 36kDa	Conjugate: -
Calculated MW: 36kDa	Isotype: IgG

Immunogen Information

Background: Enables transcription coactivator activity. Involved in histone H3 acetylation and histone deubiquitination. Located in nucleoplasm. Part of SAGA complex and transcription factor TFIIIC complex.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: SPT3; SPT3L; SUPT3H

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 50-130 of human SUPT3H (NP_852001.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.