

CAB18271

Product Information

Product SKU:	CAB18271	Gene ID:	116039	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	38kDa	Conjugate:	-
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background: OSR2 is a mammalian homolog of the Drosophila odd-skipped family of transcription factors (Lan et al., 2004 [PubMed 15175245]).

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: OSR2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-312 of human OSR2 (NP_001135934.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.