

CAB18443

Product Information

Product SKU:	CAB18443	Gene ID:	25980	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	43kDa	Conjugate:	-
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes the homolog of the yeast A1-alpha2 repressin protein that is involved in mRNA splicing. Alternately spliced transcript variants have been found for this gene.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: CGI-23; C20orf4; AAR2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human AAR2 (NP_056326.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.