

GPR168 Rabbit Polyclonal Antibody



CAB18532

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

38KDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

116535

Uniprot

Q96AM1

Synonyms:

GPR140; GPR168; MRGF; RTA; MRGPRF

Immunogen:

Recombinant protein of human GPR168.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - GPR168 Rabbit pAb (CAB18532)

