

## CAB18781

---

### Product Information

<b>Product SKU:</b>	CAB18781	<b>Gene ID:</b>	5546	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	66kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	52kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a protein that may play a role in pre-mRNA splicing. Chromosomal translocations (X;1)(p11;q21) that result in fusion of this gene to TFE3 (GeneID 7030) have been associated with papillary renal cell carcinoma. A PRCC-TFE3 fusion protein is expressed in affected carcinomas and is likely associated with altered gene transactivation. This fusion protein has also been associated with disruption of the cell cycle.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** TPRC; RCCP1; Papillary renal cell carcinoma

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide of Human Papillary renal cell carcinoma.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.