

CAB18800

Product Information

Product SKU:	CAB18800	Gene ID:	5315	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	58kDa	Conjugate:	-
Calculated MW:	58kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein involved in glycolysis. The encoded protein is a pyruvate kinase that catalyzes the transfer of a phosphoryl group from phosphoenolpyruvate to ADP, generating ATP and pyruvate. This protein has been shown to interact with thyroid hormone and may mediate cellular metabolic effects induced by thyroid hormones. This protein has been found to bind Opa protein, a bacterial outer membrane protein involved in gonococcal adherence to and invasion of human cells, suggesting a role of this protein in bacterial pathogenesis. Several alternatively spliced transcript variants encoding a few distinct isoforms have been reported.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: PK3; TCB; p58; OIP3; PKM2; CTHBP; THBP1; HEL-S-30; PKM1-specific

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human PKM1-specific (NP_001193728.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.