

CAB18804

Product Information

Product SKU:	CAB18804	Gene ID:	76421	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	20kDa	Conjugate:	-
Calculated MW:	21kDa	Isotype:	IgG

Immunogen Information

Background: Involved in reciprocal meiotic recombination and synaptonemal complex assembly. Located in chromosome. Is expressed in ovary; spermatid; spermatocyte; spermatogonium; and testis. Orthologous to human C1orf146 (chromosome 1 open reading frame 146).

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: SCRE; Spo16; C1orf146

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 47-173 of mouse C1orf146 (NP_877422.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.