

# 5-Methylcytosine (5mC) Rabbit Polyclonal Antibody



**CAB18805**

---

## Product Information

<b>Product SKU:</b> CAB18805	<b>Gene ID:</b> -	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Species independent

---

## Additional Information

<b>Observed MW:</b> -	<b>Conjugate:</b> -
<b>Calculated MW:</b> -	<b>Isotype:</b> IgG

---

## Immunogen Information

**Background:**

**Recommended Dilution:** DB, 1:500 - 1:2000

**Synonyms:** -

**Purification Method:** Affinity purification

**Immunogen:** 5mC

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Contact Details | Dublin, Ireland**

**Email:** [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | **Web:** [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.