

CAB18835

Product Information

Product SKU:	CAB18835	Gene ID:	64760	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	82kDa	Conjugate:	-
Calculated MW:	82kDa	Isotype:	IgG

Immunogen Information

Background: Able to activate MAPK/ERK and TGFB signaling pathways. May regulate the activity of genes involved in intestinal barrier function and immunoprotective inflammation (By similarity. May play a role in cell proliferation.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200 IP,0.5µg-4µg antibody for 400µg-600µg extracts of whole cells

Synonyms: FAI16; RAI16; FAM160B2; RAM160B2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 74-303 of human FAM160B2 (NP_073586.5).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.