

CAB18836

Product Information

Product SKU:	CAB18836	Gene ID:	239170	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	82kDa	Conjugate:	-
Calculated MW:	82kDa	Isotype:	IgG

Immunogen Information

Background: Orthologous to human FHIP2B (FHF complex subunit HOOK interacting protein 2B).

Recommended Dilution: WB, 1:500 - 1:1000 IHC-P, 1:50 - 1:200 IF/ICC, 1:50 - 1:200 IP, 0.5µg-4µg antibody for 400µg-600µg extracts of whole cells

Synonyms: Rai16; Fam160b2; G430067P06Rik; FAM160B2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 74-303 of mouse FAM160B2 (NP_919326.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.