

CAB18856

Product Information

Product SKU:	CAB18856	Gene ID:	9410	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	40kDa	Conjugate:	-
Calculated MW:	39kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a component of the U5 small nuclear ribonucleoprotein (snRNP) particle. The U5 snRNP is part of the spliceosome, a multiprotein complex that catalyzes the removal of introns from pre-messenger RNAs.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: 40K; SPF38; WDR57; PRP8BP; HPRP8BP; PRPF8BP; SNRNP40

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-357 of human SNRNP40 (NP_004805.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.