

CAB19151

Product Information

Product SKU:	CAB19151	Gene ID:	4340018	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	66kDa	Conjugate:	-
Calculated MW:	66kDa	Isotype:	IgG

Immunogen Information

Background: Required for the synthesis of amylose in endosperm.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: Wx; WX1; Wx-b; waxy; GBSS1; WX-mq; WX1-1; Wx-oz; GBSS-I; OsGBSS1; 134P10.7; OsGBSS I

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsGBSS I .

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.